



# ozOptics

shop.ozoptics.com  
www.ozoptics.com

219 Westbrook Road  
Ottawa, ON, Canada, K0A 1L0

Toll free: 1-800-361-5415  
Telephone: 1-613-831-0981  
Fax: 1-613-836-5089  
sales@ozoptics.com

## IN-LINE MULTIMODE FIBER SPECKLE HOMOGENIZER

### Features


- Over 90% speckle contrast reduction
- In-Line fiber device- zero optical loss
- Plug-and-play integration
- Supports a wide range of wavelengths and optical power levels
- Low power consumption
- USB type C power input
- Compact enclosure- OEM Integration and Laboratory use



### Applications


- Medical endoscopic and optical coherence imaging
- Flow cytometry and DNA sequencing
- Interferometry and Fluorescence Microscopy
- Bioanalytical instrumentation
- Machine vision and laser projection
- Speckle reduction in fiber-delivered laser systems

More than 90% speckle reduction



### Product Description


This high-performance optical Speckle Homogenizer is designed to significantly reduce speckle noise in multimode fiber laser systems. Its compact, plug-and-play design integrates both the optics and driving electronics within a single enclosure, simplifying installation and system integration. By delivering high-efficiency speckle suppression in an all-in-one in-line format, it enables superior imaging quality and system reliability across a broad spectrum of applications. With unmatched versatility, it is the go-to solution for optical engineers and system designers looking to unlock the full potential of multimode fiber imaging technology.



## Optical and Electrical Specifications



Part number: FHOM-11-IRVIS-200/240-QM-33-3A-1	
Parameter	Typical Value
Wavelength	400 nm – 1550 nm
Fiber core size	100 µm – 400 µm
Fiber core type	Round, Square, Hexagonal
Power Input	Standard USB port C 5V DC, 0.5A
Dimensions	7.6 x 3.1 x 3.1 cm
Weight	130 g
Customizable fiber and connector types See Oz Standard tables <a href="https://www.ozoptics.com/ALLNEW_PDF/DTS0079.pdf">https://www.ozoptics.com/ALLNEW_PDF/DTS0079.pdf</a> or Contact Us	

## Improvement in Speckle Contrast



$$C(\text{Speckle Contrast}) = \frac{\sigma(\text{standard dev of intensity})}{I(\text{average intensity})}$$

## 200/240 µm Round Core Fiber at 635 nm Wavelength and 100Hz Frequency



## Speckle Homogenizer FHOM-11-W-a/b-F-XY-JD-L

**W:** For multimode fibers specify either IRVIS for visible and infrared applications (400–2000 nm), or UVVIS for ultraviolet and visible applications (200–700 nm)

**a/b:** = Fiber core and cladding diameters, in microns:  
See tables 1 to 5 of the Standard Tables for standard fiber sizes.  
[https://www.ozoptics.com/ALLNEW\\_PDF/DTS0079.pdf](https://www.ozoptics.com/ALLNEW_PDF/DTS0079.pdf)

**F:** Fiber Type:  
M = Multimode fiber  
QM = High power multimode fiber

**L:** Patchcord length, in meters

**JD** = Jacket Diameter:

3 = 3 mm OD PVC loose tube with Kevlar  
3A = 3 mm OD armored  
3AS = 3 mm OD stainless steel armored  
5A = 5 mm OD armored  
5AS = 5 mm OD stainless steel armored

See table 7 of the Standard Tables for drawings  
[https://www.ozoptics.com/ALLNEW\\_PDF/DTS0079.pdf](https://www.ozoptics.com/ALLNEW_PDF/DTS0079.pdf)

**XY:** = Input and Output Connector Types:  
Refer to Table 6 of the Standard Tables Data Sheet DTS0079.  
[https://www.ozoptics.com/ALLNEW\\_PDF/DTS0079.pdf](https://www.ozoptics.com/ALLNEW_PDF/DTS0079.pdf)

**Example: FHOM-11-IRVIS-200/240-QM-33-3A-1**